

Product Name: LAYN Mouse Monoclonal Antibody**Catalog #: AMM83593**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 43 kDa

Antigen Information


Gene Name	LAYN
Alternative Names	LAYN
Gene ID	
SwissProt ID	Q6UX15
Immunogen	Purified recombinant fragment of human LAYN

Background**Research Area**


Image Data



Lane 1: 15 μ g of MCF7 lysate; Lane 2: 15 μ g of U-251MG lysate; Lane 3: 15 μ g of K-562 lysate; Result: LAYN Mouse Monoclonal Antibody can detect LAYN by Western blotting;



Lane 1: K-562 lysate; Lane 2: LAYN immunoprecipitated from K-562 lysate by AMM83593; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83593 can immunoprecipitate LAYN;



IHC-P analysis of stomach tissue by LAYN antibody (AMM83593)