
Product Name: S100A12 Mouse Monoclonal Antibody**Catalog #: AMM83595**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200
Molecular Weight	10.6kDa

Antigen Information


Gene Name	S100A12
Alternative Names	S100A12
Gene ID	
SwissProt ID	P80511
Immunogen	Purified recombinant fragment of human S100A12

Background**Research Area**

Image Data



Result: S100A12 Mouse Monoclonal Antibody can detect human S100A12 by Western blotting.



IHC-P analysis of spleen tissue by S100A12 antibody (AMM83595)