
Product Name: CFB Mouse Monoclonal Antibody**Catalog #: AMM83598**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	86kDa

Antigen Information


Gene Name	CFB
Alternative Names	CFB
Gene ID	
SwissProt ID	P00751
Immunogen	Purified recombinant fragment of human CFB

Background**Research Area**

Image Data



Western blot analysis using human CFB mouse mAb against human plasma



Lane 1: CFB immunoprecipitated from plasma by AMM83598; Lane 2: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83598 can immunoprecipitate CFB.