

Product Name: LMOD1 Mouse Monoclonal Antibody
Catalog #: AMM83605



Summary

Production Name	LMOD1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	LMOD1
Alternative Names	LMOD1
Gene ID	
SwissProt ID	P29536.Purified recombinant fragment of human LMOD1

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200,IHC:1:100-1:200
Molecular Weight	67kDa


Background

Product Name: LMOD1 Mouse Monoclonal Antibody
Catalog #: AMM83605




Research Area


Image Data



Western blot analysis using LMOD1 mouse mAb against rat uterus tissue lysate



Lane 1: Rat uterus tissue lysate; Lane 2: LMOD1 immunoprecipitated from rat uterus tissue lysate by AMM83605; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83605 can immunoprecipitate LMOD1;



IHC-P analysis of endometrial smooth muscle tissue by LMOD1 antibody (AMM83605)

Product Name: LMOD1 Mouse Monoclonal Antibody
Catalog #: AMM83605



Note

For research use only.