
Product Name: ASRGL1 Mouse Monoclonal Antibody**Catalog #: AMM83608**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	ASRGL1
Alternative Names	ASRGL1
Gene ID	
SwissProt ID	Q7L266
Immunogen	Purified recombinant fragment of human ASRGL1

Background**Research Area**

Image Data



Lane 1: Hep G2 lysate; Lane 2: ASRGL1 immunoprecipitated from Hep G2 lysate by ASRGL1 Mouse Monoclonal Antibody; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: ASRGL1 Mouse Monoclonal Antibody can immunoprecipitate ASRGL1;