
Product Name: VAMP1 Mouse Monoclonal Antibody**Catalog #: AMM83612**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200
Molecular Weight	12.9kDa

Antigen Information


Gene Name	VAMP1
Alternative Names	VAMP1
Gene ID	
SwissProt ID	P23763
Immunogen	Purified recombinant fragment of human VAMP1

Background**Research Area**

Image Data



Western blot analysis using VAMP1 mouse mAb against rat brain tissue lysate



IHC-P analysis of cerebral cortex tissue by anti-VAMP1 antibody (AMM83612)