
Product Name: SYNPR Mouse Monoclonal Antibody**Catalog #: AMM83615**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	SYNPR
Alternative Names	SYNPR
Gene ID	
SwissProt ID	Q8TBG9
Immunogen	Purified recombinant fragment of human SYNPR

Background**Research Area**

Image Data



IHC-P analysis of cerebral cortex tissue by SYNPR antibody (SYNPR Mouse Monoclonal Antibody)