
Product Name: DMXL2 Mouse Monoclonal Antibody**Catalog #: AMM83617**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	339.6kDa

Antigen Information


Gene Name	DMXL2
Alternative Names	DMXL2
Gene ID	
SwissProt ID	Q8TDJ6
Immunogen	Purified recombinant fragment of human DMXL2

Background**Research Area**

Image Data



Lane 1: 15 μ g of Jurkat lysate; Lane 2: 15 μ g of HeLa lysate; Lane 3: 15 μ g of HEK-293 lysate; Result: DMXL2 Mouse Monoclonal Antibody can detect DMXL2 by Western blotting;



Lane 1: HEK-293 lysate; Lane 2: DMXL2 immunoprecipitated from HEK-293 lysate by AMM83617; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83617 can immunoprecipitate DMXL2;