
Product Name: DMXL2 Mouse Monoclonal Antibody**Catalog #: AMM83618**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	DMXL2
Alternative Names	DMXL2
Gene ID	
SwissProt ID	Q8TDJ6
Immunogen	Purified recombinant fragment of human DMXL2

Background**Research Area**

Image Data



Lane 1: HEK-293 lysate; Lane 2: DMXL2 immunoprecipitated from HEK-293 lysate by DMXL2 Mouse Monoclonal Antibody; Lane 3: The same Lane 2 but use IgG isotype control antibody; Result: DMXL2 Mouse Monoclonal Antibody can immunoprecipitate DMXL2;



IHC-P analysis of cerebral cortex tissue by DMXL2 antibody (AMM83618)