
Product Name: TSNAX Mouse Monoclonal Antibody**Catalog #: AMM83622**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:250
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 33 kDa

Antigen Information


Gene Name	TSNAX
Alternative Names	TSNAX
Gene ID	
SwissProt ID	Q99598
Immunogen	Purified recombinant fragment of human TSNAX

Background**Research Area**

Image Data



Lane 1: 15 µg of MCF7 lysate; Lane 2: 15 µg of HeLa lysate; Lane 3: 15 µg of HEK-293 lysate; Lane 4: 15 µg of SK-OV-3 lysate; Result: TSNAX Mouse Monoclonal Antibody can detect TSNAX by Western blotting.



IHC-P analysis of testis tissue by TSNAX antibody (AMM83622)