

Product Name: ASPM Mouse Monoclonal Antibody
Catalog #: AMM83623



Summary

Production Name	ASPM Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ASPM
Alternative Names	ASPM
Gene ID	
SwissProt ID	Q8IZT6.Purified recombinant fragment of human ASPM

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:250
Molecular Weight	382kDa


Background

Product Name: ASPM Mouse Monoclonal Antibody
Catalog #: AMM83623




Research Area

Image Data



Western blot analysis using ASPM mouse mAb against HEK-293 lysate



IHC-P analysis of adrenal gland tissue by ASPM antibody (AMM83623)

Note

For research use only.