

Product Name: JCHAIN Mouse Monoclonal Antibody**Catalog #: AMM83624**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:200-1:1000,IHC 1:200-1:500
Molecular Weight	18.1kDa

Antigen Information


Gene Name	JCHAIN
Alternative Names	JCHAIN
Gene ID	
SwissProt ID	P01591
Immunogen	Purified recombinant fragment of human JCHAIN

Background**Research Area**

Image Data



Western blot analysis using human JCHAIN mouse mAb against human plasma lysate



IHC-P analysis of tonsil tissue by JCHAIN antibody (P154999)