
Product Name: SLC32A1 Mouse Monoclonal Antibody**Catalog #: AMM83628**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:200
Molecular Weight	57.4kDa

Antigen Information


Gene Name	SLC32A1
Alternative Names	SLC32A1
Gene ID	
SwissProt ID	Q9H598
Immunogen	Purified recombinant fragment of human SLC32A1

Background**Research Area**

Image Data



Western blot analysis using SLC32A1 mouse mAb against rat brain tissue



IHC-P analysis of cerebral cortex tissue by anti-SLC32A1 antibody (AMM83628)