

Product Name: GMPR2 Mouse Monoclonal Antibody**Catalog #: AMM83629**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	37.9kDa

Antigen Information


Gene Name	GMPR2
Alternative Names	GMPR2
Gene ID	
SwissProt ID	Q9P2T1
Immunogen	Purified recombinant fragment of human GMPR2

Background**Research Area**

Image Data



Western blot analysis using GMPR2 mouse mAb against Hep G2 lysate



Lane 1: GMPR2 immunoprecipitated from Hep G2 lysate by AMM83629; Lane 2: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83629 can immunoprecipitate GMPR2.