

Product Name: PSMD1 Mouse Monoclonal Antibody**Catalog #: AMM83630**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:200-1:500
Molecular Weight	105.8kDa

Antigen Information


Gene Name	PSMD1
Alternative Names	PSMD1
Gene ID	
SwissProt ID	Q99460
Immunogen	Purified recombinant fragment of human PSMD1

Background**Research Area**

Image Data



Western blot analysis using PSMD1 mouse mAb against Hep G2 lysate



IHC-P analysis of testis tissue by PSMD1 antibody (AMM83630)