
Product Name: RGS18 Mouse Monoclonal Antibody**Catalog #: AMM83632**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:50-1:200
Molecular Weight	/

Antigen Information

Gene Name	RGS18
Alternative Names	RGS18
Gene ID	
SwissProt ID	Q9NS28
Immunogen	Purified recombinant fragment of human RGS18

Background**Research Area**

Image Data



IHC-P analysis of kidney tissue by RGS18 antibody (RGS18 Mouse Monoclonal Antibody)