

Product Name: ARAF Mouse Monoclonal Antibody
Catalog #: AMM83634



Summary

Production Name	ARAF Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ARAF
Alternative Names	ARAF
Gene ID	
SwissProt ID	P10398. Purified recombinant fragment of human ARAF

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	68kDa


Background

Product Name: ARAF Mouse Monoclonal Antibody
Catalog #: AMM83634



Research Area

Image Data



Lane 1: 15 μ g of K-562 lysate Lane 2: 15 μ g of HeLa lysate Lane 3: 15 μ g of HEK-293 lysate

Note

For research use only.