

Product Name: PGF Mouse Monoclonal Antibody
Catalog #: AMM83640



Summary

Production Name	PGF Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	PGF
Alternative Names	PGF
Gene ID	
SwissProt ID	P49763. Purified recombinant fragment of human PGF

Application

Dilution Ratio	IHC:1:50-1:200
Molecular Weight	/


Background

Product Name: PGF Mouse Monoclonal Antibody
Catalog #: AMM83640



Research Area

Image Data



IHC-P analysis of placenta tissue by PGF antibody (AMM83640)

Note

For research use only.