

Summary

Production Name	FETUB Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	FETUB
Alternative Names	FETUB
Gene ID	
SwissProt ID	Q9UGM5.Purified recombinant fragment of human FETUB

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:50-1:200
Molecular Weight	42kDa


Background

Product Name: FETUB Mouse Monoclonal Antibody
Catalog #: AMM83642




Research Area

Image Data



Western blot analysis using FETUB mouse mAb against kidney tissue lysate



IHC-P analysis of hepatocellular carcinoma tissue by FETUB antibody (AMM83642)

Note

For research use only.