

Product Name: TRIM21 Mouse Monoclonal Antibody**Catalog #: AMM83645**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:200
Molecular Weight	Calculated MW: 54 kDa; Observed MW: 50 kDa

Antigen Information


Gene Name	TRIM21
Alternative Names	TRIM21
Gene ID	
SwissProt ID	P19474
Immunogen	Purified recombinant fragment of human TRIM21

Background**Research Area**

Image Data



Western blot analysis using TRIM21 mouse mAb against HL-60 lysate



IHC-P analysis of tonsil tissue by TRIM21 antibody (AMM83645)