
Product Name: DNAJB1 Mouse Monoclonal Antibody**Catalog #: AMM83654**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:200,IP 1:100-1:200
Molecular Weight	45kDa

Antigen Information


Gene Name	DNAJB1
Alternative Names	DNAJB1
Gene ID	
SwissProt ID	P25685
Immunogen	Purified recombinant fragment of human DNAJB1

Background**Research Area**


Image Data



Western blot analysis using DNAJB1 mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate; Lane 2: DNAJB1 immunoprecipitated from Jurkat lysate by AMM83654; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83654 can immunoprecipitate DNAJB1;



IHC-P analysis of esophagus tissue by DNAJB1 antibody (AMM83654)