
Product Name: STAT2 Mouse Monoclonal Antibody**Catalog #: AMM83661**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	113kDa

Antigen Information


Gene Name	STAT2
Alternative Names	STAT2
Gene ID	
SwissProt ID	P52630
Immunogen	Purified recombinant fragment of human STAT2

Background**Research Area**

Image Data



Western blot analysis using STAT2 mouse mAb against K-562 lysate



Lane 1: K-562 lysate; Lane 2: STAT2 immunoprecipitated from K-562 lysate by AMM83661; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83661 can immunoprecipitate STAT2;