

**Product Name: PLCB2 Mouse Monoclonal Antibody**  
**Catalog #: AMM83662**



---

## Summary

|                        |                                 |
|------------------------|---------------------------------|
| <b>Production Name</b> | PLCB2 Mouse Monoclonal Antibody |
| <b>Description</b>     | Mouse Monoclonal Antibody       |
| <b>Host</b>            | Mouse                           |
| <b>Application</b>     | WB                              |
| <b>Reactivity</b>      | Human                           |

## Performance

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| <b>Conjugation</b>  | Unconjugated                                                                             |
| <b>Modification</b> | Unmodified                                                                               |
| <b>Isotype</b>      | Mouse IgG1                                                                               |
| <b>Clonality</b>    | Monoclonal                                                                               |
| <b>Form</b>         | Liquid                                                                                   |
| <b>Storage</b>      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| <b>Buffer</b>       | Purified antibody in PBS with 0.05% sodium azide                                         |
| <b>Purification</b> | Affinity Purification                                                                    |

## Immunogen

|                          |                                                      |
|--------------------------|------------------------------------------------------|
| <b>Gene Name</b>         | PLCB2                                                |
| <b>Alternative Names</b> | PLCB2                                                |
| <b>Gene ID</b>           |                                                      |
| <b>SwissProt ID</b>      | Q00722. Purified recombinant fragment of human PLCB2 |

## Application

|                         |                 |
|-------------------------|-----------------|
| <b>Dilution Ratio</b>   | WB:1:500-1:1000 |
| <b>Molecular Weight</b> | 134kDa          |

## Background


**Product Name: PLCB2 Mouse Monoclonal Antibody**  
**Catalog #: AMM83662**



---

## Research Area

## Image Data



Western blot analysis using PLCB2 mouse mAb against THP-1 lysate

## Note

For research use only.