

Product Name: IRAG2 Mouse Monoclonal Antibody**Catalog #: AMM83665**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:20-1:100,IP 1:100-1:200
Molecular Weight	Calculated MW::Observed MW:61kD

Antigen Information


Gene Name	IRAG2
Alternative Names	IRAG2
Gene ID	
SwissProt ID	Q12912
Immunogen	Purified recombinant fragment of human IRAG2

Background**Research Area**


Image Data



Western blot analysis using IRAG2 mouse mAb against Ramos lysate



Lane 1: Ramos lysate ; Lane 2: IRAG2 immunoprecipitated from Ramos lysate by AMM83665 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83665 can immunoprecipitate IRAG2;



IHC-P analysis of tonsil tissue by IRAG2 antibody (AMM83665)