

Product Name: OLR1 Mouse Monoclonal Antibody**Catalog #: AMM83668**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:20-1:100
Molecular Weight	/

Antigen Information

Gene Name	OLR1
Alternative Names	OLR1
Gene ID	
SwissProt ID	P78380
Immunogen	Purified recombinant fragment of human OLR1

Background**Research Area**

Image Data



IHC-P analysis of kidney tissue by OLR1 antibody (OLR1 Mouse Monoclonal Antibody)