

Summary

Production Name	MYH1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	ICC
Reactivity	Human,Mouse,Rat,Fruitfly,Nematode

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%BSA and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	MYH1
Alternative Names	fast skeletal muscle; MYH3; HEMHC; Myosin-1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1
Gene ID	4619.0
SwissProt ID	P12882. Synthetic Peptide of Myosin Heavy Chain

Application

Dilution Ratio	ICC:1:50-1:200
Molecular Weight	/

Background


Muscle contraction.

Product Name: MYH1 Mouse Monoclonal Antibody
Catalog #: AMM84935



Research Area

Image Data



Immunofluorescence analysis of MYH1 in Immunofluorescence Staining of nematode tissue using MYH1 antibody.

Note

For research use only.