

Product Name: Rabbit IgG Fc Mouse Monoclonal Antibody
Catalog #: AMM84972

Summary

Production Name	Rabbit IgG Fc Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Rabbit

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Rabbit IgG Fc
Alternative Names	Antibody; Ig gamma chain C region; Ig kappa-b4 chain C region; Ig kappa-b5 chain C region; Ig kappa-b9 chain C region; Ig lambda chain C region; K-BAS.
Gene ID	-
SwissProt ID	-Recombinant Protein of Rabbit IgG (Fc specific, Heavy Chain)

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 55 kDa; Observed MW: 55 kDa


Product Name: Rabbit IgG Fc Mouse Monoclonal Antibody
Catalog #: AMM84972

Background

Constant region of immunoglobulin heavy chains. Immunoglobulins, also known as antibodies, are membrane-bound or secreted glycoproteins produced by B lymphocytes.

Research Area

Image Data



Western blot analysis of Reduced Rabbit IgG, H+L using Rabbit IgG Fc antibody.

Note

For research use only.