

Product Name: Prealbumin Mouse Monoclonal Antibody**Catalog #: AMM84983**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa

Antigen Information


Gene Name	Prealbumin
Alternative Names	Prealbumin; TBPA; Amyloidosis I; Transthyretin
Gene ID	7276.0
SwissProt ID	P02766
Immunogen	Purified recombinant protein expressed in E.coli.

Background


Thyroid hormone-binding protein. Probably transports thyroxine from the bloodstream to the brain.

Research Area


Image Data



Western blot analysis of Prealbumin in Hela, Jurkat, mouse brain, MOLT4, Human serum lysates using Prealbumin antibody



Immunohistochemistry analysis of paraffin-embedded Human Liver Carcinoma Tissue using Prealbumin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using Prealbumin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.