

Product Name: DDX3 Mouse Monoclonal Antibody
Catalog #: AMM85097



Summary

Production Name	DDX3 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB, ICC, IP
Reactivity	Human, Mouse, Rat, Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	DDX3
Alternative Names	DBX; DDX3; HLP2; DDX14.
Gene ID	1654.0
SwissProt ID	O00571. Purified recombinant human DDX3 protein fragments expressed in E.coli.

Application

Dilution Ratio	WB 1:500-1:1000, ICC 1:50-1:200, IP 1:20-1:50
Molecular Weight	Calculated MW: 73 kDa; Observed MW: 73 kDa

Background

Product Name: DDX3 Mouse Monoclonal Antibody
Catalog #: AMM85097




Multifunctional ATP-dependent RNA helicase.

Research Area

Apoptosis

Image Data



Western blot analysis of DDX3 in COS7, Jurkat, C6, 3T3 and HeLa lysates using DDX3 antibody.

Note

For research use only.