

Summary

Production Name	MIPEP Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	MIPEP
Alternative Names	Mitochondrial intermediate peptidase, MIP, MIPEP, MIP
Gene ID	4285.0
SwissProt ID	Q99797. This MIPEP monoclonal antibody is generated from mouse immunized with MIPEP recombinant protein.

Application

Dilution Ratio	WB:1:100-1:1000
Molecular Weight	80.4kDa

Background


Cleaves proteins, imported into the mitochondrion, to their mature size.

Product Name: MIPEP Mouse Monoclonal Antibody
Catalog #: AMM85914



Research Area

Image Data



MIPEP Antibody (Cat. #AMM85914) western blot analysis in 293 cell line lysates (35µg/lane). This demonstrates the MIPEP antibody detected the MIPEP protein (arrow).

Note

For research use only.