

Product Name: UBL5(Ascites) Mouse Monoclonal Antibody**Catalog #: AMM85930**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:8000
Molecular Weight	8.5kDa

Antigen Information


Gene Name	UBL5(Ascites)
Alternative Names	Ubiquitin-like protein 5, UBL5
Gene ID	59286.0
SwissProt ID	Q9BZL1
Immunogen	Purified His-tagged UBL5 protein(Fragment) was used to produced this monoclonal antibody.

Background

Ubiquitin-like proteins (UBLs) are thought to be reversible modulators of protein function rather than protein degraders like ubiquitin (MIM 191339).

Research Area

Image Data



UBL5 Antibody western blot analysis in mouse heart tissue lysates (35µg/lane). This demonstrates the UBL5 antibody detected the UBL5 protein (arrow).