
Product Name: TPIPb(Ascites) Mouse Monoclonal Antibody**Catalog #: AMM85946**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:2000-1:5000
Molecular Weight	61.1kDa


Antigen Information

Gene Name	TPIPb(Ascites)
Alternative Names	Phosphatidylinositol 3,5-trisphosphate 3-phosphatase TPTE2, Lipid phosphatase TPIP, TPTE and PTEN homologous inositol lipid phosphatase, TPTE2, TPIP
Gene ID	93492.0
SwissProt ID	Q6XPS3
Immunogen	Purified His-tagged TPIPb protein was used to produced this monoclonal antibody.

Background

Research Area

Image Data



TPIPb Antibody (C-term) western blot analysis in 293,A431,A2058,mouse NIH/3T3 cell line lysates (35 μ g/lane).This demonstrates the TPIPb antibody detected the TPIPb protein (arrow).