
Product Name: GMFG(Ascites) Mouse Monoclonal Antibody**Catalog #: AMM85949**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:10000-1:40000
Molecular Weight	16.8kDa


Antigen Information

Gene Name	GMFG(Ascites)
Alternative Names	Glia maturation factor gamma, GMF-gamma, GMFG
Gene ID	9535.0
SwissProt ID	O60234
Immunogen	Purified His-tagged GMFG protein was used to produced this monoclonal antibody.


Background**Research Area**

TGF-beta signaling pathway

Image Data



GMFG Antibody (N-term) western blot analysis in Jurkat cell line lysates (35µg/lane). This demonstrates the GMFG antibody detected the GMFG protein (arrow).



GMFG Antibody (Cat. #AMM85949) western blot analysis in Jurkat and MDA-MB-453 cell line lysates (35µg/lane). This demonstrates the GMFG antibody detected the GMFG protein (arrow).