

Summary

Production Name	EEF1E1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,ICC,FC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	EEF1E1
Alternative Names	Eukaryotic translation elongation factor 1 epsilon-1, Aminoacyl tRNA synthetase complex-interacting multifunctional protein 3, Elongation factor p18, Multisynthase complex auxiliary component p18, EEF1E1, AIMP3, P18
Gene ID	9521.0
SwissProt ID	O43324. This EEF1E1 antibody is generated from a mouse immunized with a recombinant protein of human EEF1E1.

Application

Dilution Ratio	WB:1:2002-1:4000,ICC:1:25,FC:1:25
Molecular Weight	19.8kDa

Product Name: EEF1E1 Mouse Monoclonal Antibody
Catalog #: AMM85989




Background

Positive modulator of ATM response to DNA damage.

Research Area

Image Data



All lanes : Anti-EEF1E1 Antibody at 1:2000-1:4000 dilution

Note

For research use only.