
Product Name: DPYSL5 Mouse Monoclonal Antibody**Catalog #: AMM85990**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ICC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:25-1:50
Molecular Weight	61.4kDa

Antigen Information


Gene Name	DPYSL5 Dihydropyrimidinase-related protein 5, DRP-5, CRMP3-associated molecule, CRAM,
Alternative Names	Collapsin response mediator protein 5, CRMP-5, UNC33-like phosphoprotein 6, ULIP-6, DPYSL5, CRMP5, ULIP6
Gene ID	56896.0
SwissProt ID	Q9BPU6
Immunogen	This DPYL5 antibody is generated from a mouse immunized with a recombinant protein of human DPYL5.

Background

May have a function in neuronal differentiation and/or axon growth.

Research Area

Image Data



All lanes : Anti-DPYL5 Antibody at 1:4000 dilution