

Summary

Production Name	ALDOC Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG3
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	ALDOC
Alternative Names	Fructose-bisphosphate aldolase C, Brain-type aldolase, ALDOC, ALDC
Gene ID	230.0
SwissProt ID	P09972. Purified His-tagged ALDOC protein was used to produced this monoclonal antibody.

Application

Dilution Ratio	WB:1:2000,FC:1:25
Molecular Weight	39.5kDa


Background

Product Name: ALDOC Mouse Monoclonal Antibody
Catalog #: AMM85995



Research Area

Image Data



All lanes : Anti-ALDOC Antibody (C-term) at 1:2000 dilution

Note

For research use only.