
Product Name: KIF1B Mouse Monoclonal Antibody**Catalog #: AMM85996**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	130kDa

Antigen Information

Gene Name	KIF1B
Alternative Names	KIF1B
Gene ID	
SwissProt ID	O60333-3
Immunogen	This KIF1B antibody is generated from a mouse immunized with a recombinant protein.

Background**Research Area**

Image Data



All lanes : Anti-KIF1B Antibody at 1:2000 dilution