

Product Name: RAB3A Mouse Monoclonal Antibody
Catalog #: AMM85998



Summary

Production Name	RAB3A Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB, ICC, FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	RAB3A
Alternative Names	Ras-related protein Rab-3A, RAB3A
Gene ID	5864.0
SwissProt ID	P20336. This RAB3A antibody is generated from a mouse immunized with a recombinant protein.

Application

Dilution Ratio	WB 1:500-1:2000, ICC 1:25-1:50, FC 1:25-1:50
Molecular Weight	30.0kDa

Background


Product Name: RAB3A Mouse Monoclonal Antibody
Catalog #: AMM85998



Involved in exocytosis by regulating a late step in synaptic vesicle fusion. Could play a role in neurotransmitter release by regulating membrane flow in the nerve terminal.

Research Area

Image Data



All lanes : Anti-RAB3A Antibody at 1:2000 dilution

Note

For research use only.