
Product Name: RAB3B Mouse Monoclonal Antibody**Catalog #: AMM86000**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,FC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:500,FC 1:25-1:50
Molecular Weight	24.8kDa

Antigen Information


Gene Name	RAB3B
Alternative Names	Ras-related protein Rab-3B, RAB3B
Gene ID	5865.0
SwissProt ID	P20337
Immunogen	This RAB3B antibody is generated from a mouse immunized with a recombinant protein of human RAB3B.

Background

Protein transport. Probably involved in vesicular traffic (By similarity).

Research Area

Image Data



All lanes : Anti-RAB3B Antibody at 1:2000 dilution