
Product Name: CDX1 Mouse Monoclonal Antibody**Catalog #: AMM86005**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	28.1kDa

Antigen Information

Gene Name	CDX1
Alternative Names	Homeobox protein CDX-1, Caudal-type homeobox protein 1, CDX1
Gene ID	1044.0
SwissProt ID	P47902
Immunogen	Purified His-tagged CDX1 protein was used to produced this monoclonal antibody.

Background

Could play a role in the terminal differentiation of the intestine.

Research Area

Image Data

All lanes : Anti-CDX1 Antibody (C-term) at 1:1000-1:2000 dilution

