

Product Name: MKRN2 Mouse Monoclonal Antibody
Catalog #: AMM86006



Summary

Production Name	MKRN2 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	MKRN2
Alternative Names	Probable E3 ubiquitin-protein ligase makorin-2, 632-, RING finger protein 62, MKRN2, RNF62
Gene ID	23609.0
SwissProt ID	Q9H000.This MKRN2 antibody is generated from a mouse immunized with recombinant protein.

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:500,ICC 1:25-1:50,FC 1:25-1:50
Molecular Weight	46.9kDa

Product Name: MKRN2 Mouse Monoclonal Antibody
Catalog #: AMM86006




Background

E3 ubiquitin ligase catalyzing the covalent attachment of ubiquitin moieties onto substrate proteins.

Research Area

Image Data



All lanes : Anti-MKRN2 Antibody at 1:4000 dilution

Note

For research use only.