
Product Name: RBM6 Mouse Monoclonal Antibody**Catalog #: AMM86007**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	128.6kDa

Antigen Information


Gene Name	RBM6
Alternative Names	RNA-binding protein 6, Lung cancer antigen NY-LU-12, Protein G16, RNA-binding motif protein 6, RNA-binding protein DEF-3, RBM6, DEF3
Gene ID	10180.0
SwissProt ID	P78332
Immunogen	This RBM6 antibody is generated from a mouse immunized with a recombinant protein of human RBM6.

Background

Specifically binds poly(G) RNA homopolymers in vitro.

Research Area

Image Data



All lanes : Anti-RBM6 Antibody at 1:1000 dilution