
Product Name: RAB27A Mouse Monoclonal Antibody**Catalog #: AMM86013**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,FC
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000,FC 1:25-1:50
Molecular Weight	24.9kDa

Antigen Information


Gene Name	RAB27A
Alternative Names	Ras-related protein Rab-27A, Rab-27, GTP-binding protein Ram, RAB27A, RAB27
Gene ID	5873.0
SwissProt ID	P51159
Immunogen	This RAB27A antibody is generated from a mouse immunized with a recombinant protein of human RAB27A.

Background

Plays a role in cytotoxic granule exocytosis in lymphocytes. Required for both granule maturation and granule docking and priming at the immunologic synapse.

Research Area

Image Data



Anti-RAB27A Antibody at 1:2000 dilution + K562 whole cell lysate