
Product Name: RRAS2 Mouse Monoclonal Antibody**Catalog #: AMM86014**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	24.0kDa

Antigen Information

Gene Name	RRAS2
Alternative Names	Ras-related protein R-Ras2, Ras-like protein TC21, Teratocarcinoma oncogene, RRAS2, TC21
Gene ID	22800.0
SwissProt ID	P62070
Immunogen	This RRAS2 antibody is generated from a mouse immunized with a recombinant protein of human RRAS2.


Background

It is a plasma membrane-associated GTP-binding protein with GTPase activity. Might transduce growth inhibitory signals across the cell membrane, exerting its effect through an effector shared with the Ras proteins but in an antagonistic fashion.

Research Area

TGF-beta signaling pathway, MAPK signaling pathway

Image Data



All lanes : Anti-RRAS2 Antibody at 1:4000 dilution