
Product Name: RAB5B Mouse Monoclonal Antibody**Catalog #: AMM86016**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	23.7kDa

Antigen Information


Gene Name	RAB5B
Alternative Names	Ras-related protein Rab-5B, RAB5B
Gene ID	5869.0
SwissProt ID	P61020
Immunogen	This RAB5B antibody is generated from a mouse immunized with a recombinant protein of human RAB5B.

Background

Protein transport. Probably involved in vesicular traffic (By similarity).

Research Area

Image Data



Anti-RAB5B Antibody at 1:2000 dilution + A431 whole cell lysate