

Product Name: PBRM1 Mouse Monoclonal Antibody**Catalog #: AMM86028**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	192.9kDa

Antigen Information

Gene Name	PBRM1
Alternative Names	Protein polybromo-1, hPB1, BRG1-associated factor 180, BAF180, Polybromo-1D, PBRM1, BAF180, PB1
Gene ID	55193.0
SwissProt ID	Q86U86
Immunogen	This PBRM1 antibody is generated from a mouse immunized with a recombinant protein of human PBRM1.


Background

Involved in transcriptional activation and repression of select genes by chromatin remodeling (alteration of DNA-nucleosome

topology). Acts as a negative regulator of cell proliferation.

Research Area

Image Data



Anti-PBRM1 Antibody at 1:1000 dilution + HT-1080 whole cell lysate