

Product Name: RAB1B Mouse Monoclonal Antibody**Catalog #: AMM86041**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	22.2kDa

Antigen Information

Gene Name	RAB1B
Alternative Names	Ras-related protein Rab-1B, RAB1B
Gene ID	81876.0
SwissProt ID	Q9H0U4
Immunogen	This RAB1B antibody is generated from a mouse immunized with a recombinant protein between 1-201 amino acids from human RAB1B.

Background


The small GTPases Rab are key regulators of intracellular membrane trafficking, from the formation of transport vesicles to their fusion with membranes. Rabs cycle between an inactive GDP-bound form and an active GTP-bound form that is able to recruit

to membranes different set of downstream effectors directly responsible for vesicle formation, movement, tethering and fusion. RAB1B regulates vesicular transport between the endoplasmic reticulum and successive Golgi compartments. Plays a role in the initial events of the autophagic vacuole development which take place at specialized regions of the endoplasmic reticulum.

Research Area

Autophagy

Image Data



All lanes : Anti-RAB1B Antibody at 1:4000 dilution