

Product Name: FGFR1OP2 Mouse Monoclonal Antibody**Catalog #: AMM86044**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	29.4kDa

Antigen Information


Gene Name	FGFR1OP2
Alternative Names	FGFR1 oncogene partner 2, FGFR1OP2
Gene ID	26127.0
SwissProt ID	Q9NVK5
Immunogen	This FGFR1OP2 antibody is generated from a mouse immunized with a recombinant protein between 1-172 amino acids of human FGFR1OP2.

Background

May be involved in wound healing pathway.

Research Area

Image Data



All lanes : Anti-FGFR1OP2 Antibody at 1:4000 dilution