

Summary

Production Name	ELAVL1 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB, ICC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	ELAVL1
Alternative Names	ELAV-like protein 1, Hu-antigen R, HuR, ELAVL1, HUR
Gene ID	1994.0
SwissProt ID	Q15717. This ELAVL1 monoclonal antibody is generated from mouse immunized with ELAVL1 recombinant protein.

Application

Dilution Ratio	WB 1:500-1:2000, ICC 1:20-1:50
Molecular Weight	36.1kDa

Background


Product Name: ELAVL1 Mouse Monoclonal Antibody
Catalog #: AMM86049



The protein encoded by this gene is a member of the ELAVL protein family. This encoded protein contains 3 RNA-binding domains and binds cis-acting AU-rich elements. It destabilizes mRNAs and thereby regulates gene expression.

Research Area

Image Data



All lanes : Anti- at 1:1000 dilution

Note

For research use only.