

Product Name: RAB20 Mouse Monoclonal Antibody
Catalog #: AMM86053



Summary

Production Name	RAB20 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	RAB20
Alternative Names	Ras-related protein Rab-20, RAB20
Gene ID	55647.0
SwissProt ID	Q9NX57. This RAB20 antibody is generated from a mouse immunized with recombinant protein from human RAB20.

Application

Dilution Ratio	WB:1:2000-1:4000
Molecular Weight	26.3kDa

Background

Plays a role in apical endocytosis/recycling. Plays a role in the maturation and acidification of phagosomes that engulf


Product Name: RAB20 Mouse Monoclonal Antibody
Catalog #: AMM86053



pathogens, such as S.aureus and M.tuberculosis. Plays a role in the fusion of phagosomes with lysosomes.

Research Area

Image Data



All lanes : Anti-RAB20 Antibody at 1:2000-1:4000 dilution

Note

For research use only.